

Oral hygiene

Dental Prosthesis Cleaning System C2S PRO®

The C2S Pro cleaning system was developed to improve dental care in healthcare facilities. It significantly lowers bacteria and fungi (candidiasis) levels, reducing the occurrence of some systemic diseases (upper respiratory disease). Regular care of dental prostheses improves the general health of residents. In addition, daily use of the C2S Pro reduces referrals to specialists (dentists, physicians and prosthodontists).

The C2S Pro cleaning system, developed for daily and restorative cleaning of all types of removable and screw-retained dentures through dynamic brushing with special alloy bristles, helps clean dentures thoroughly and prevents the formation of plaque and tartar.



Hygienic



Ergonomic

Key points

- ▶ Ease of use and makes cleaning easier
- ▶ Eliminate plaque that naturally lines dentures
- ▶ Simple and effective for ensuring the long life of dentures
- ▶ Limits the use of specialists (dentists, dental technicians, etc.)

Oral hygiene

Dental Prosthesis Cleaning System C2S Pro®

Cleaning system

- ▶ Dim: L. 20 x W. 17,5 x H. 19,5 cm
- ▶ Cleaning operation: magnets rotating at 1800 rpm
- ▶ Timing: 2 min / 5 min / 10 min / 15 min / 20 min
- ▶ Automatic system shutdown
- ▶ Certification: CE
- ▶ Power supply: separate CE-certified transformer, input : 100-240 V / output : 12 V
- ▶ Warranty: 2 years



Description	Reference
Dental Prosthesis Cleaning System C2S PRO®	2CPRO

Removable battery

- ▶ Type: Lithium Polymère
- ▶ Dimensions: 19,4 x 12 x 2,1 cm
- ▶ Colour: silver grey
- ▶ Indicators: blue LED
- ▶ Autonomy: 7 h
- ▶ INPUT power supply: 19 V / 1,2 A
- ▶ OUTPUT power supply:
 - DC 12 V / 2,5 A (for C2S system)
 - DC 16 V / 4 A
 - DC 19 V / 4 A
 - USB 5 V / 1,2 A
- ▶ Capacity: 30000 mAh



Description	Reference
Removable battery C2S BAT®	2CPROBAT

Cleaning box

- ▶ Dim : Ø 11 x H. 9,5 cm
- ▶ Certification: CE
- ▶ Reporting: name labels



Description	Reference
2 cleaning pink boxes	2CBOLRM

Bicarbonate tablets

- ▶ Effervescent tablets for daily use
- ▶ Composition :
 - Sodium bicarbonate
 - Cider vinegar
 - Thyme-mint fragrance
 - Magnesium stearate
 - Microcrystalline cellulose
 - Blue colouring (E133)



Description	Reference
Box of 200 tablets C2S TAB®	2CTAB

Alloy bristles

Description	Reference
Box of 20 blisters alloy bristles	2CPINS

