

Neop | 400 x 400 double-column distribution trolleys

Daily distribution of medicines with pill dispensers and nominative trays

- Robustness: excellent resistance, one-piece structure made of aluminium
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Customization thanks to a wide range of drawers and accessories
- Ward name identification
- Ergonomic and handy: designed for medical staff comfort
- Unique design and harmonious colours



Robust



Hygienic



Ergonomic



Customizable

- Double-column trolley
- Height: 100 cm
- Roller shutter door or door wing
- Code lock (with spare key) or badge lock
- 2 available inside panel models (slides) according to internal fittings
- Ergonomic gripping area

Customizable structure

3 available colours

Magenta Pink

Apple Green

Sky Blue



Trolley configuration

STEP 1 - Trolley structure

Choice of trolley structure according to inside panel models

Trolley Height 100 cm*	
Standard inside panels 40 mm level	Inside panels for telescopic drawers
8Y244H100P4	8Y244H100T

This reference includes:

- **Structure**
- **Wheels**
- **Inside panels**

STEP 2 - Trolley Customization

Selection of elements 1 to 4 is compulsory

Closure system, locking system and colour strips must be double ordered for a double-column trolley

			Trolley with standard inside panels ref. 8Y244H100P4	Trolley with inside panels for telescopic drawers ref. 8Y244H100T
1	Upper shelf	compact laminate	8Y001_4400CP2	8Y001_4400CP2
		Roller shutter door	8Y002_4010F3	-
2	Closing system	Door wing for telescopic drawers	-	8Y002_0010F2
		For roller shutter door	8Y003CODE3	-
3	Code lock	For door wing for telescopic drawers	-	8Y003CODE2
		For roller shutter door	8Y003BADGE3	-
	Badge lock	For door wing for telescopic drawers	-	8Y003BADGE2
		Magenta Pink	8Y004_0010RM	8Y004_0010RM
4	Colour	Apple Green	8Y004_0010VP	8Y004_0010VP
		Sky Blue	8Y004_0010BC	8Y004_0010BC

*For precise height of working surface, see technical data



Compact laminate upper shelf



Roller shutter closure



Code lock

Magenta Pink

Apple Green



Door wing



Badge lock

Sky Blue

STEP 3 - Accessories and internal fittings

Extension shelf, standard rail, ergonomic handle, motorization, arm for laptop, upper frame, drawer...

Please refer to product data sheets.



Technical data

Overall dimensions of the trolley



Useful dimensions of the worktop



Worktop height

H. 1064 mm

Materials

Structure	Complete welded one-piece aluminium structure, EPOXY anti-corrosion treated
Work surface	Compact laminate, 10 mm thickness
Slides	ABS (Acrylonitrile Butadiene Styrene) OR ballbearing slide made of galvanized steel with shock absorber
Closing system	Roller shutter door: PP (Polypropylene), 2 mm thickness Door wing: Aluminium with EPOXY painting
Wheels	Simple non marking rubber wheels - Polyamide (PA) cover
Bumpers	4 non marking round stock blocks made of EVA (Ethylene Vinyl Acetate)
Strip of colour	70 micron calendered PVC vinyl film with acrylic adhesive

Additional information

Wheels	125 mm diameter - 2 breaking wheels + 2 swivelling wheels Wheel and wheel-pin double blocking
Empty weight	40 kg
Capacity	- 38 levels of 40 mm for trolley with roller shutter door - 40 levels of 40 mm for trolley without roller shutter door - 10 telescopic drawers H 80 mm
Locking	Code lock / Badge lock

Packaging

Packing	Individually on pallet strapped and film-wrapped or on wheels with proper protection	Dimensions L X D x H	80 x 120 x 120 cm
Gross weight (trolley + packaging)	52 kg	Volume	1,15 m ³