

Neop | 400 x 400 single-column nursing trolleys

Care management system

- Robustness: excellent resistance, one-piece structure made of aluminium
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Customization thanks to a wide range of drawers and accessories
- Ergonomic and handy: designed for medical staff comfort
- Ward name identification



Robust



Hygienic



Ergonomic



Customizable

- Single-column trolley
- Height* 100 cm
- Roller shutter door or door wing
- Code lock or badge lock
- 2 available inside panel models (slides) according to required internal fittings

Customizable structure

3 available colours

Magenta pink

Apple green

Sky blue



*For precise height of working surface, see « Technical data »

Standard configurations - PREMIUM

Structure WITH internal fittings and peripheral accessories.
Remember to add the color reference of the trolley.

TELESCOPIC DRAWERS WITHOUT CLOSING SYSTEM


3 small / 2 medium / 1 large

400 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- 3 telescopic drawers H. 80 mm for ISO trays H. 50 mm
- 2 telescopic drawers H. 160 mm for ISO trays H. 100 mm
- 1 telescopic drawer H. 240 mm for ISO tray H. 100 mm

Reference

8YS44H100TT4



TELESCOPIC DRAWERS WITH DOOR WING

3 small / 2 medium / 1 large

400 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 3 telescopic drawers H. 80 mm for ISO trays H. 50 mm
- 2 telescopic drawers H. 160 mm for ISO trays H. 100 mm
- 1 telescopic drawer H. 240 mm for ISO tray H. 100 mm
- 3 standard rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 1 double bag holder
- 1 large side basket

Reference

8YS44H100TT3



TELESCOPIC DRAWERS WITH DOOR WING

5 medium

400 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 5 telescopic drawers H. 160 mm for ISO trays H. 100 mm
- 3 standard rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 1 double bag holder
- 1 large side basket

Reference

8YS44H100TT1



TELESCOPIC DRAWERS WITH DOOR WING


1 small / 3 medium / 1 large

400 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Door wing
- Code lock
- 1 telescopic drawer H. 80 mm for ISO trays H. 50 mm
- 3 telescopic drawers H. 160 mm for ISO trays H. 100 mm
- 1 telescopic drawer H. 240 mm for ISO trays H. 100 mm
- 3 standard rails (top right, bottom right, bottom left)
- 1 standard support for sharp container
- 1 double bag holder
- 1 large side basket

Reference

8YS44H100TT2



JOINT DRAWERS WITH ROLLER SHUTTER DOOR

1 petit / 3 moyens / 1 grand

400 x 400 nursing trolley H100 equipped with:

- Thermoplastic upper shelf
- Roller shutter door
- Code lock
- 1 drawer 400 x 400 H. 75 mm
- 3 drawers 400 x 400 H. 115 mm
- 1 drawer 400 x 400 H. 195 mm
- 3 standard rails (top left, bottom left, bottom right)
- 1 foldable extension shelf

Reference

8YS44H100TJ1



Standard configurations – BASIC

Structure WITH internal fittings and peripheral accessories.

Remember to add the color reference of the trolley.


TROLLEY WITHOUT FOLDABLE EXTENSION SHELF

400 x 400 nursing trolley H100 equipped with:

- Inside panels for telescopic drawers
- Thermoplastic upper shelf
- Door wing
- Code lock
- 3 standard rails (top right, bottom right, bottom left)

Reference

8YCHST144H100




TROLLEY WITHOUT FOLDABLE EXTENSION SHELF

400 x 400 nursing trolley H100 equipped with:

- Inside panels for telescopic drawers
- Thermoplastic upper shelf
- Door wing
- Code lock
- 2 standard rails (bottom right, middle left)
- 1 standard support for sharp container
- 1 simple bag holder

Reference

8YCHST2144H100




TROLLEY WITH FOLDABLE EXTENSION SHELF

400 x 400 nursing trolley H100 equipped with:

- Inside panels for telescopic drawers
- Thermoplastic upper shelf
- Door wing
- Code lock
- 2 standard rails (bottom right, bottom left)
- 1 foldable extension shelf

Reference

8YCHST144H100T



Custom configuration

STEP 1 - Trolley structure

Choice of trolley structure according to inside panel models

Trolley height 100 cm*	
Standard inside panels 40 mm levels	Inside panels for telescopic drawers
8Y144H100P4	8Y144H100T

This reference includes:

- **Structure**
- **Wheels**
- **Inside panels**

STEP 2 - Trolley customization

Selection of elements 1 to 4 is compulsory

Trolley with standard
inside panels
ref. 8Y144H100P4

Trolley with inside panels
for telescopic drawers
ref. 8Y144H100T

1	Upper shelf	thermoplastic	8Y001_4400TH	8Y001_4400TH
2	Closing system	Roller shutter door	8Y002_4010F3	-
		Door wing for telescopic drawers	-	8Y002_0010F2
3	Code lock	For roller shutter door	8Y003CODE3	-
		For door wing for telescopic drawers	-	8Y003CODE2
	Badge lock	For roller shutter door	8Y003BADGE3	-
		For door wing for telescopic drawers	-	8Y003BADGE2
4	Colour	Magenta pink	8Y004_0010RM	8Y004_0010RM
		Apple green	8Y004_0010VP	8Y004_0010VP
		Sky blue	8Y004_0010BC	8Y004_0010BC

*For precise height of working surface, see next page



Accessories and internal fittings

Extension shelf, standard rail, ergonomic handle, arm for laptop, upper frame, drawer...

Please refer to product data sheets.



Technical data

Overall dimensions of the trolley



Useful dimensions of the worktop



Thermoplastic upper shelf
with ergonomic gripping area

Worktop height

H. 1057 mm

Materials

Structure	Complete welded one-piece aluminium structure, EPOXY anti-corrosion treated
Work surface	Thermoplastic: ABS / PMMA
Slides	ABS (Acrylonitrile Butadiene Styrene) or ballbearing slide made of galvanized steel with shock absorber
Closing system	Roller shutter door: PP (Polypropylene), 2 mm thickness Door wing: Aluminium with EPOXY painting
Wheels	Simple non marking rubber wheels - Polyamide cover (PA)
Bumpers	4 non marking round stop blocks made of EVA (Ethylene Vinyl Acetate)
Strip of colour	70 micron calendered PVC vinyl film, with acrylic adhesive

Additional information

Wheels	125 mm diameter - 2 breaking wheels + 2 swivelling wheels Wheel and wheel-pin double blocking
Empty weight	20 kg
Capacity	- 19 levels of 40 mm for trolley with roller shutter door - 20 levels of 40 mm for trolley without roller shutter door - 10 telescopic drawers H 80 mm
Locking system	Code lock (with spare key) / Badge lock

Packaging

Packing	Individually on pallet strapped and film-wrapped, or on wheels, with proper protection	Dimensions L X D x H	80 x 120 x 120 cm
Gross weight (trolley + packaging)	32 kg	Volume	1,15 m ³