

### Storage of medicines and medical devices

- Robustness: excellent resistance, one-piece structure made of steel
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Customization thanks to a wide range of drawers and accessories
- Ward name identification
- Ergonomic and handy: designed for medical staff comfort
- Unique design and harmonious colours



Robust



Hygienic



Ergonomic



Customizable

- Single-column cabinet
- 2 heights: 160 or 180 cm
- Roller shutter door, or without closure
- Code lock (with spare key) or badge lock, or without lock
- Height 180 cm: horizontal or tilted levels in upper position

#### Customizable structure

3 available colours

Magenta pink

Apple green

Sky blue



## Standard configurations - PREMIUM

Structure WITH internal fittings


### MD – ISO TRAY

400 x 600 single-column mobile storage cabinet H180 equipped with:

- Horizontal levels on entire height
- Roller shutter door
- Code lock
- 6 perforated ISO trays H100
- 3 perforated ISO trays H200
- 9 pairs of end stops for ISO tray

Reference	Description
8YA146H180DMBAL_RM	NEOP "Premium" 1 col. 400x600 mobile storage cabinet H180 - MD – ISO tray - magenta pink
8YA146H180DMBAL_BC	NEOP "Premium" 1 col. 400x600 mobile storage cabinet H180 - MD - ISO tray - sky blue
8YA146H180DMBAL_VP	NEOP "Premium" 1 col. 400x600 mobile storage cabinet H180 - MD - ISO tray - apple green

*Dividers, label holders and label covers to order if necessary*



## Standard configurations – BASIC

Structure WITHOUT internal fittings


### Mobile storage cabinets with tilted levels in upper position

#### MOBILE CABINET WITH ROLLER SHUTTER DOOR, WITHOUT LOCK

400 x 600 single-column mobile storage cabinet H180 equipped with:

- Tilted levels in upper position
- Roller shutter door

Reference	Description
8YAH146XH180FR_RM	NEOP "Basic" 1 col. 400x600 mobile storage cabinet H180 Tilted levels without lock - magenta pink
8YAH146XH180FR_BC	NEOP "Basic" 1 col. 400x600 mobile storage cabinet H180 Tilted levels without lock - sky blue
8YAH146XH180FR_VP	NEOP "Basic" 1 col. 400x600 mobile storage cabinet H180 Tilted levels without lock - apple green




#### MOBILE CABINET WITH ROLLER SHUTTER DOOR AND CODE LOCK

400 x 600 single-column mobile storage cabinet H180 equipped with:

- Tilted levels in upper position
- Roller shutter door
- Code lock

Reference	Description
8YAH146XH180FRC_RM	NEOP "Basic" 1 col. 400x600 mobile storage cabinet H180 Tilted levels + code lock - magenta pink
8YAH146XH180FRC_BC	NEOP "Basic" 1 col. 400x600 mobile storage cabinet H180 Tilted levels + code lock - sky blue
8YAH146XH180FRC_VP	NEOP "Basic" 1 col. 400x600 mobile storage cabinet H180 Tilted levels + code lock - apple green



## Mobile cabinets with horizontal levels on entire height

### MOBILE CABINET WITHOUT ROLLER SHUTTER DOOR

400 x 600 single-column mobile storage cabinet H180 equipped with horizontal levels on entire height

Reference	Description
8YAH146H180_RM	NEOP "Basic" mobile storage cabinet 1 col. 400x600 H180 without roller shutter door - magenta pink
8YAH146H180_BC	NEOP "Basic" mobile storage cabinet 1 col. 400x600 H180 without roller shutter door - sky blue
8YAH146H180_VP	NEOP "Basic" mobile storage cabinet 1 col. 400x600 H180 without roller shutter door - apple green



### MOBILE CABINET WITH ROLLER SHUTTER DOOR WITHOUT LOCK

400 x 600 single-column mobile storage cabinet H180 equipped with:

- Horizontal levels on entire height
- Roller shutter door

Reference	Description
8YAH146H180FR_RM	NEOP "Basic" mobile storage cabinet 1 col. 400x600 H180 roller shutter door without lock - magenta pink
8YAH146H180FR_BC	NEOP "Basic" mobile storage cabinet 1 col. 400x600 H180 roller shutter door without lock - sky blue
8YAH146H180FR_VP	NEOP "Basic" mobile storage cabinet 1 col. 400x600 H180 roller shutter door without lock - apple green




### MOBILE CABINET WITH ROLLER SHUTTER DOOR AND CODE LOCK

400 x 600 single-column mobile storage cabinet H180 equipped with:

- Horizontal levels on entire height
- Roller shutter door
- Code lock

Reference	Description
8YAH146H180FRC_RM	NEOP "Basic" mobile storage cabinet 1 col. 400x600 H180 roller shutter door + code lock - magenta pink
8YAH146H180FRC_BC	NEOP "Basic" mobile storage cabinet 1 col. 400x600 H180 roller shutter door + code lock - sky blue
8YAH146H180FRC_VP	NEOP "Basic" mobile storage cabinet 1 col. 400x600 H180 roller shutter door + code lock - apple green



## Custom configurations

### STEP 1 - Cabinet structure

Choice of cabinet structure according to inside panel models

	Cabinet height 160 cm* with bumpers	Cabinet height 180 cm* with bumpers	Cabinet height 180 cm* without bumpers
Horizontal levels on entire height 40 mm level	8Y146H160P4	8Y146H180P4	8Y146H180SPCP4
Tilted levels in upper position 40 mm level	-	8Y146H180PI	8Y146H180SPCPI

This reference includes:

- **Structure**
- **Wheels**
- **Inside panels**

### STEP 2 - Cabinet customization

Selection of elements 1 to 3

	Cabinet height 160 cm*	Cabinet height 180 cm*
<b>1</b> <b>Roller shutter closure</b>	8Y002_4616F3	8Y002_4618F3

Storage cabinet can be ordered without closure.

**BE CAREFUL !** For a cabinet with horizontal levels on entire height without roller shutter door, you have to add a end stop ref. [8Y014](#).

<b>2</b>	<b>Code lock</b>	8Y003CODE3	8Y003CODE3
	<b>Badge lock</b>	8Y003BADGE3	8Y003BADGE3

Storage cabinet can be ordered without lock.

<b>3</b>	<b>Colour</b>	Magenta pink: 8Y004_0016RM	Magenta pink: 8Y004_0018RM
		Apple green: 8Y004_0016VP	Apple green: 8Y004_0018VP
		Sky blue: 8Y004_0016BC	Sky blue: 8Y004_0018BC

\*For precise height, see technical data



Roller shutter closure



Code lock



Badge lock

Magenta pink

Apple green

Sky blue


## Accessories and internal fittings

Standard rail, push handle, trays...

Please refer to product data sheets



## Technical data



### Materials

<b>Structure</b>	Complete welded one-piece structure, electro-galvanized steel, EPOXY anti-corrosion treated
<b>Slides</b>	ABS (Acrylonitrile Butadiene Styrene)
<b>Closing system</b>	Roller shutter door: PP (Polypropylene), 2 mm thickness
<b>Wheels</b>	Simple non marking rubber wheels - Polyamide cover (PA)
<b>Bumpers</b>	4 non marking round stop blocks made of EVA (Ethylene Vinyl Acetate)
<b>Strip of colour</b>	70 micron calendered PVC vinyl film, with acrylic adhesive

### Additional information

<b>Wheels</b>	125 mm diameter - 2 breaking wheels + 2 swivelling wheels wheel and wheel-pin double blocking	
<b>Empty weight</b>	H 160: 54 kg	H 180: 58 kg
<b>Capacity</b>	<b>H 160:</b> - 35 levels of 40 mm for cabinet with roller shutter door - 36 levels of 40 mm for cabinet without roller shutter door	<b>H 180 with horizontal levels:</b> - 39 levels of 40 mm for cabinet with roller shutter door - 40 levels of 40 mm for cabinet without roller shutter door  <b>H 180 with tilted levels :</b> 7 tilted 40 mm levels + 25 horizontal 40 mm levels
<b>Locking system</b>	Code lock / Badge lock	

### Packaging

<b>Packing</b>	Individually on pallet strapped and film-wrapped, or on wheels, with proper protection	<b>Dimensions L X D x H</b>	H 160: 80 x 120 x 185 cm H 180: 80 x 120 x 200 cm
<b>Gross weight (cabinet + packaging)</b>	H 160: 66 kg H 180: 70 kg	<b>Volume</b>	1,90 m <sup>3</sup>