

Neop | 600 x 400 emergency trolleys with joint drawers

Material and products immediately available for first aid

- Robustness: excellent resistance, one-piece structure made of aluminium
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Customization thanks to a wide range of drawers and accessories
- Ward name identification
- Ergonomic and handy: designed for medical staff comfort
- Unique design and harmonious colours



Robust



Hygienic



Ergonomic



Customizable

- Height 100 cm
- Ergonomic gripping area
- Door wing lock with 1 batch of 20 seals. Access to the content of the trolley breaking the seal by simple twist. After each use, the trolley is reequipped, then closed by a single-use numbered seal.
- 5 or 6 joint drawers
- Cardiac massage board on support

Colour

Garance red



Standard configuration - PREMIUM


EMERGENCY TROLLEYS – 5 joint drawers

Peripheral accessories

- 1 batch of 20 seals
- 1 cardiac massage board holder
- 1 cardiac massage board
- 1 monitoring holder
- 1 IV pole with 2 hooks
- 1 lateral standard rail - bottom right
- 1 oxygen tank holder

Internal fittings

- 1 drawer 600 x 400 height 75 mm
- 3 drawers 600 x 400 height 155 mm
- 1 drawer 600 x 400 height 235 mm



Reference	Description
8Y164H100P4.U2	600 x 400 emergency trolleys - 5 joint drawers


EMERGENCY TROLLEYS – 6 joint drawers

Peripheral accessories

- 1 batch of 20 seals
- 1 cardiac massage board holder
- 1 cardiac massage board
- 1 monitoring holder
- 1 IV pole with 2 hooks
- 1 lateral standard rail - bottom right
- 1 oxygen tank holder

Internal fittings

- 3 drawers 600 x 400 height 75 mm
- 2 drawers 600 x 400 height 155 mm
- 1 drawer 600 x 400 height 235 mm



Reference	Description
8Y164H100P4.U3	600 x 400 emergency trolleys - 6 joint drawers

Accessories and internal fittings



Ergonomic handle, upper frame, standard rail, bag holder...
Please refer to product data sheets.

Technical data

Overall dimensions of the trolley



Les cotes hors tout sont relevées avec les roues dans la position la moins encombrante.

Useful dimensions of the worktop



Worktop height

H. 1057 mm

Materials

Structure	Complete welded and EPOXY anti-corrosion treated one-piece structure
Work surface	ABS / PMMA thermoplastic
Closing system	Door wing: Aluminium with EPOXY painting
Wheels	Simple non marking rubber wheels - Polyamide cover (PA)
Bumpers	4 non marking round stop blocks made of EVA (Ethylene Vinyl Acetate)
Strip of colour	70 micron calendered PVC vinyl film with acrylic adhesive

Additional information

Wheels	125 mm diameter - 2 breaking wheels + 2 swivelling wheels - wheel and wheel-pin double blocking
Empty weight	24 kg
Capacity	20 levels of 40 mm
Locking system	Seals

Packaging

Packing	Individually on pallet strapped and film-wrapped, or on wheels, with proper protection	Dimensions L X D x H	L. 80 x D. 120 x H. 145 cm
Gross Weight (trolley + packaging)	36 kg	Volume	1,15 m ³