

Neop | 600 x 400 transfer cabinets

Supply of medicines or medical devices

- Robustness: excellent resistance, torsion-resistant steel frame, primer powder + polyester treated
- Towing for the transport is possible on the dedicated locations on the 4 corners, or anywhere on the external transfer frame
- 100 % hygiene: designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Customization thanks to a wide range of drawers and accessories
- Ergonomic and handy: designed for medical staff comfort
- Unique design and harmonious colours



Robust



Hygienic



Ergonomic



Customizable

- Single-column transfer cabinet
- 2 heights: 160 or 180 cm
- Roller shutter door
- Code lock (with spare key) or badge lock
- 2 available inside panel models (slides) according to required internal fittings

Customizable structure

3 available colours

Magenta pink

Apple green

Sky blue



Standard configurations

Transfer cabinets H160 with standard inside panels - 40 mm levels

Reference	Description
8YATR164H160_RM_V2	NEOP "Basic" transfer cabinet 1 col. 600x400 H160 roller shutter door without lock - magenta pink
8YATR164H160_VP_V2	NEOP "Basic" transfer cabinet 1 col. 600x400 H160 roller shutter door without lock - apple green
8YATR164H160_BC_V2	NEOP "Basic" transfer cabinet 1 col. 600x400 H160 roller shutter door without lock - sky blue

Transfer cabinets H160 with optimised inside panels - 62 mm levels

Reference	Description
8YATH160FRC_RM_ED2	NEOP Basic transfer cabinet 1 col. 600x400 H160 62 levels - roller shutter door without lock - magenta pink
8YATH160FRC_BC_ED2	NEOP Basic transfer cabinet 1 col. 600x400 H160 62 levels - roller shutter door without lock - sky blue
8YATH160FRC_VP_ED2	NEOP Basic transfer cabinet 1 col. 600x400 H160 62 levels - roller shutter door without lock - apple green

Transfer cabinets H180 with optimised inside panels - 62 mm levels

Reference	Description
8YATH180FRC_RM_ED2	NEOP Basic transfer cabinet 1 col. 600x400 H180 62 levels - roller shutter door without lock - magenta pink
8YATH180FRC_BC_ED2	NEOP Basic transfer cabinet 1 col. 600x400 H180 62 levels - roller shutter door without lock - sky blue
8YATH180FRC_VP_ED2	NEOP Basic transfer cabinet 1 col. 600x400 H180 62 levels - roller shutter door without lock - apple green

Custom configurations

STEP 1 - Cabinet structure

Choice of cabinet structure according to inside panel models

	Cabinet height 160 cm*	Cabinet height 180 cm*
Standard inside panels 40 mm level	8Y164H160P4.TR_V2	8Y164H180P4.TR_V2
Optimized inside panels* 50 mm level	8Y164H160P5.TR_V2	8Y164H180P5.TR_V2

This reference includes:

- **Structure**
- **Framework**
- **Wheels**
- **Inside panels**

* for Modulo XL pill dispensers , 600 x 400 Modulo Classic and Eurojour pill dispensers

STEP 2 - Cabinet Customization

Selection of elements 1 to 3 is compulsory

		Cabinet height 160 cm*	Cabinet height 180 cm*	
1	Roller shutter door	8Y002_6016F3	8Y002_6018F3	
2	Code lock	8Y003CODE3	8Y003CODE3	
	Badge lock	8Y003BADGE3	8Y003BADGE3	
3	Colour	Magenta pink	8Y004_0016RM	8Y004_0018RM
		Apple green	8Y004_0016VP	8Y004_0018VP
		Sky blue	8Y004_0016BC	8Y004_0018BC
4	Wheels	8Y031	8Y031	

*for precise height, see next page



Roller shutter closure



Code lock



Badge lock

Magenta pink

Apple green

Sky blue

STEP 3 - Options

Option 4 pivoting wheels including 1 directional and 2 brakes	8Y043
--	-------

Accessories and internal fittings

Please refer to product data sheets



Technical data



Materials

Structure	Complete welded one-piece structure, electro-galvanized steel, EPOXY anti-corrosion treated
Framework	Machine welded steel framework, powder primer coated + polyester treated
Slides	ABS (Acrylonitrile Butadiene Styrene)
Closing system	Roller shutter door: PP (Polypropylene), 2 mm thickness
Wheels	Semi-elastic rubber wheels - pressed steel bracket
Bumpers	4 rubber corners

Additional information

Wheels	160 mm diameter - 2 fixed wheels + 2 swivelling wheels with brake	
Empty weight (framework + structure)	Cabinet H160: 29 + 48 = 77 kg	Cabinet H180: 30 + 52 = 82 kg
Capacity	Cabinet H160 with roller shutter door: - 35 levels of 40 mm - 27 levels of 50 mm - 22 levels of 62 mm	Cabinet H180 with roller shutter door: - 39 levels of 40 mm - 30 levels of 50 mm - 24 levels of 62 mm
Locking system	Code lock / Badge lock	
Strip of colour	70 micron calendered PVC vinyl film, with acrylic adhesive	

Packaging

Packing	Individually on pallet strapped and film-wrapped, or on wheels, with proper protection	Dimensions	H 160: L. 88 x D. 80 x H. 200 cm H 180: L. 88 x D. 80 x H. 215 cm
Gross weight (cabinet + packaging)	H 160: 89 kg H 180: 94 kg	Volume	H 160: 1,4 m ³ H 180: 1,5 m ³