

Care furniture

Agily® 600 x 400 trolleys

The Agily® dressing trolleys are adaptable and can be equipped for different uses. PRATICIMA will help you to adapt your equipment according to your needs.

Thanks to their innovative design and their ultra-light components, these carts are very good tools for hospital work.

Photos are not contractual



Advantages

- ▶ Ergonomic, handy and light aluminium structure
- ▶ Adaptable according to your needs: shelves and 600 x 400 trays...
- ▶ 2 types of shelves: thermoplastic and compact
- ▶ Resistant to the disinfecting products of hospital environments
- ▶ Custom height on request




Hygienic



Ergonomic

Agily® 600 x 400 trolleys


Trolleys with compact shelf



Configuration n°11

600 x 400 trolleys equipped with a one top compact shelf with handle

> Ref. AGIL6_11




Configuration n°12

600 x 400 trolleys equipped with:

- ▶ 1 top compact shelf with handle
- ▶ 1 lower compact shelf

> Ref. AGIL6_12




Configuration n°13

600 x 400 trolleys equipped with:

- ▶ 1 top compact shelf with handle
- ▶ 1 middle compact shelf
- ▶ 1 lower compact shelf

> Ref. AGIL6_13


Trolleys with thermoplastic shelf



Configuration n°03

600 x 400 trolleys equipped with a one top thermoplastic shelf

> Ref. AGIL6_03




Configuration n°05

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 3 standard rails (top right, bottom right and top left)

> Ref. AGIL6_05




Configuration n°14

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 1 tablet holder

> Ref. AGIL6_14

Agily® 600 x 400 trolleys




Configuration n°24

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 1 drawer H. 75 mm with runners
- ▶ 3 standard rails (top right, bottom right and top left)

> Ref. AGIL6_24




Configuration n°25

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 1 drawer H. 115 mm with runners
- ▶ 4 standard rails

> Ref. AGIL6_25




Configuration n°28

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 2 drawer H. 115 mm with runners
- ▶ 4 standard rails

> Ref. AGIL6_28



Configuration n°30


600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 2 drawer H. 115 mm with runners
- ▶ 4 standard rails
- ▶ 1 plastic tray
- ▶ 1 standard support for sharp container
- ▶ 1 double bag holder

> Ref. AGIL6_30

Care furniture

Agily® 600 x 400 trolleys




Configuration n°35

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 bottle wedge 5 boxes - small model
- ▶ 1 lower thermoplastic shelf
- ▶ 1 push handle
- ▶ 1 drawer H. 75 mm with runners
- ▶ 2 drawer H. 115 mm with runners
- ▶ 1 drawer H. 155 mm with runners
- ▶ 1 compartmentalization kit 6 boxes for drawer H. 75 mm
- ▶ 1 compartmentalization kit 6 boxes for drawer H. 115 mm
- ▶ 4 standard rails
- ▶ 1 standard support for glove box
- ▶ 1 compartment sorter
- ▶ 1 standard support for sharp container
- ▶ 1 double bag holder

> Ref. AGIL6_35




Configuration n°07

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 1 standard rail (top right)
- ▶ 1 foldable extension shelf

> Ref. AGIL6_07




Configuration n°17

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 1 support pour tablette tactile numérique
- ▶ 1 standard rail (top right)
- ▶ 1 foldable extension shelf

> Ref. AGIL6_17

Agily® 600 x 400 trolleys




Configuration n°26

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 4 standard rails
- ▶ 1 foldable extension shelf
- ▶ 1 drawer H. 115 mm with runners

> Ref. AGIL6_26




Configuration n°29

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 2 drawer H. 115 mm with runners
- ▶ 4 standard rails
- ▶ 1 foldable extension shelf

> Ref. AGIL6_29




Configuration n°31

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 2 drawer H. 115 mm with runners
- ▶ 4 standard rails
- ▶ 1 foldable extension shelf
- ▶ 1 plastic tray
- ▶ 1 standard support for sharp container
- ▶ 1 double bag holder

> Ref. AGIL6_31




Configuration n°23

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 1 upper frame with 2 standard rails

> Ref. AGIL6_23

Agily® 600 x 400 trolleys




Configuration n°10

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 1 standard rail (top right)
- ▶ 1 foldable extension shelf
- ▶ 1 upper frame with 2 standard rails

> Ref. AGIL6_10




Configuration n°36

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 middle thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 3 standard rails (top right, bottom right and top left)

> Ref. AGIL6_36




Configuration n°37

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 middle thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 1 drawer H. 75 mm with runners
- ▶ 3 standard rails (top right, bottom right and top left)

> Ref. AGIL6_37



Configuration n°38

600 x 400 trolleys equipped with:

- ▶ 1 top thermoplastic shelf
- ▶ 1 middle thermoplastic shelf
- ▶ 1 lower thermoplastic shelf
- ▶ 1 drawer H. 115 mm with runners
- ▶ 2 standard rails (top right, bottom right)
- ▶ 1 standard support for sharp container
- ▶ 1 bag holder



> Ref. AGIL6_38

Care furniture

Agily® 600 x 400 trolleys

Technical data

Dimensions of trolleys with thermoplastic shelves



Care furniture

Agily® 600 x 400 trolleys

Dimensions of trolleys with compact shelves



Materials

Structure	Anodised anti-corrosive aluminium structure
Wheels	4 swivelling Ø 125 mm polyamide wheels, 1 with brake
Bumper	Non-marking rubber and buffers
Workplan	Thermoplastic : ABS / PMMA or compact laminate 10 mm thickness
Standard rail	30 x 10 mm anodized hollow aluminum profile
Drawer slides	Textured white epoxy coated 20/10 steel
Drawer	ABS PC

Various information

Wheels	4 wheels Ø 125 mm - 1 with brake and 3 swivel		
Permissible load	80 kg		
Unloaded weight	AGIL6_05 : 13,2 kg AGIL6_12 : 20,60 kg	AGIL6_13 : 25 kg AGIL6_24 : 15,4 kg	AGIL6_36 : 14,4 kg AGIL6_37 : 16,6 kg
Capacity	8 drawers 75 high		

Packaging

Conditioning	À l'unité sur palette, cerclée et filmée		
Net weight	AGIL6_05 : 23 kg AGIL6_12 : 27,4 kg	AGIL6_13 : 31,8 kg AGIL6_24 : 25 kg	AGIL6_36 : 24 kg AGIL6_37 : 26 kg
Dimensions	L. 0,88 x D. 0,80 x H. 1,10 m		
Volume	0,78 m³		













Choice of colors

<i>Pop blue</i>	<i>Apple green</i>	<i>Orange</i>	<i>Fuchsia</i>	<i>Pearl grey</i>
AGILY_BP	AGILY_VP	AGILY_OR	AGILY_FU	AGILY_GR

Care furniture

Agily® 600 x 400 trolleys

Options

Description	Picture	Reference	Description	Picture	Reference
Thermoplastic middle shelf		87194	Handle for thermoplastic Agily® trolley		87125
Thermoplastic lower shelf		87195	Pull-down extension tray (thermoplastic)		87301
Tablet intermediate in compact		87184	Wire guard on upper shelf (suction cup)		87061
Tablet lower in compact		87185	Middle wire guards on 3 sides		87068
600 x 400 trays (runners included)		TIROIR_S	Upper frame with 2 standard aluminium rails		87143
4 hooks for standard rail		28007N	Kit cale flacons 5 cases		28038