

Medicine dispensing Hebdo[®] dispensing trolleys

The Hebdo[®] range of medicine trolleys are ideal for distribution in institutions (elderly homes, nursing homes...).

Our **pill dispensers**, with fixed or removable dividers, will make storage and distribution easier for you. Their ergonomic design is an advantage for daily use.



Hygienic



Ergonomic

Key points

- ▶ Storage, transfer and distribution on one trolley only
- ▶ Securing by code lock
- ▶ Different trolley capacity
- ▶ Resistant to the disinfecting products of hospital environments

Medicine dispensing

Hebdo[®] dispensing trolleys

Trolleys

- ▶ Anodised aluminium structure
- ▶ ABS side walls, coloured body
- ▶ White thermoplastic work surface with rounded raised edges
- ▶ Curtain door with tool-free removable slats
- ▶ Code lock
- ▶ 4 polyamide swivelling wheels Ø 100 mm, 1 with brake
- ▶ Standard: 1 handle and 4 standard rails for accessories



Description	Dimensions (in mm)	Pill dispensers by trolley	Reference
Hebdo[®] 30 trolley with curtain door	L. 615 x D. 468 x H. 988	30	8D030_CODE
Hebdo[®] 60 trolley with curtain door	L. 1065 x D. 468 x H. 988	60	8D060_CODE

Compatible pill dispenses and accessories

Description	Dimensions (in mm)	Reference
1 Pill dispenser with fixed dividers	L. 319 x D. 199 x H. 42	505001BL
2 Pill dispenser with removable dividers	L. 319 x D. 199 x H. 42	50610
Batch of 10 label sheets (10 labels by sheet)		50288
3 White bin	L. 320 x D. 200 x H. 150	42016

