



Daily distribution of medicines with tracked and secure pill dispensers

- Robustness : excellent resistance, natural anodised structure, ABS-mass tinted side walls
- 100 % Hygiene : designed with a minimum of retention zones, and resistant to disinfecting products of hospital environment
- Customization thanks to a wide range of drawers and accessories.
- Ergonomic and handy : designed for medical staff comfort
- Unique design



Robust



Hygienic



Ergonomic



Customizable

Colour

Magenta Pink

- Single or double-column trolley
- Height* 100 or 120 cm
- Roller shutter door
- Code lock (with spare key)
- Standard equipment : 1 push handle and 4 standard rails for accessories

	Single-column Pili cart	Double-column Pili cart
H. 100	8Y144.PILI	8Y244.PILI
H. 120	8Y144H120.PILI	8Y244H120.PILI



Accessories and internal fittings

Extension shelf, arm for laptop, upper frame, waste bin, drawer...

Please refer to product data sheets.



* For precise height of working surface, see technical data

Technical data

Overall dimensions of the trolley



Useful dimensions of the worktop



Materials

Structure	Anodised aluminium profile
External panels	ABS (Acrylonitrile Butadiene Styrene)
Work surface	ABS / PMMA thermoplastic
Slides	ABS (Acrylonitrile Butadiene Styrene)
Closing system	Roller shutter: PP (Polypropylene)
Wheels	Simple non marking rubber wheels - Polyamide cover (PA)
Bumpers	4 non marking round stop blocks made of EVA (Ethylene Vinyl Acetate)

Additional Information

Wheels	100 mm diameter - 2 breaking wheels + 3 swivelling wheels	
Empty weight	Simple column: 21 kg / Double column : 31 kg Simple column: 23 kg / Double column : 34 kg	
Loading capacity	5 kg per tray	
Capacity	Simple column: 10 trays for PILI or Modulo Classic pill dispensers Simple column: 12 trays for PILI or Modulo Classic pill dispensers	Double column: 20 trays for PILI or Modulo Classic pill dispensers Double column: 24 trays for PILI or Modulo Classic pill dispensers
Locking system	Code / Badge	
External panel colour	Raspberry	
Equipment	Roller shutter + push handle + 4 lateral rails	

Packaging

Packing	Individually on pallet strapped and film-wrapped, or on wheels, with proper protection	Dimensions L X D x H	Simple column: 80 x 88 x 120 cm Double column: 80 x 120 x 120 cm
Gross weight (trolley + packaging)	Simple column: H.100 : 31 kg Double column: H.100 : 41 kg Simple column: H.120 : 33 kg Double column: H.120 : 44 kg	Volume	1 m ³