

Compression range

# Pinel® restraint system

The **Pinel®** design allows you to address any known restraint requirement with ONE system. At the emergency end of restraining, the aggressive patient can be secured to a seven-point restraint level in about 10 seconds ; but, adjusted safely to a comfortable one-point restraint as he/she calms. At the passive end of restraint needs, confused patients can be secured with full mobility within the bed or chair to the level that patients are not even aware of being restrained.

This single system addresses all possible restraint needs within or between (transport) any institution, e.g. Long Term Care, Emergency Room, Psychiatric, ICU or Correctional facilities...



## Advantages

- ▶ Confort: Polyester + Kevlar
- ▶ Security: magnetic locking system
- ▶ Hygiene: treatment against bacteria, easy washing & fast drying
- ▶ Liberty of movements: very flexible when the patient is secured
- ▶ Rapidity: universal system

Compression range

# Pinel® restraint system

## Complet kit

The kit contains:

- ▶ 4 limb straps
- ▶ 1 waist belt
- ▶ 2 side straps
- ▶ 1 extender
- ▶ 1 belt for trunk
- ▶ 1 pelvic strap
- ▶ 1 strap belt
- ▶ 3 keys
- ▶ 13 axis + bolt with string
- ▶ 2 axis + bolt
- ▶ 1 emergency cutter

Equipped with a color code, this case permits you to put easily the different equipments of the Pinel® range.

Description	Reference
<b>Complet kit with case</b>	SAA50100
<b>Complet kit without case</b>	SAA50200
<b>Case only</b>	SAA50900


## Bolts and accessories

Every Pinel® fixation is equipped with very effective locking system which can quickly be applied and easily used. The locking system with axis + bolt enables to block until 4 straps together. The system unlocks by putting the magnetic key in the bolt.

Resistance: 500 kg.

Description	Reference
<b>Magnetic key</b>	SAA50910
<b>Magnetic bolt</b>	SAA50920
<b>Stainless steel axis</b>	SAA50930
<b>Axis + bolt</b>	SAA50940
<b>Axis + bolt + key</b>	SAA50950

## Pelvic strap



This abdominal belt allows the patient to lateralize and move up and down. For patients who are in danger of sliding in a bed or on a chair.

The kit contains 1 abdominal belt, 2 free side straps, 1 pelvic strap, 2 axis + bolt with string, 2 axis + bolt and 3 keys.

Description	Reference
<b>Abdominal belt kit small size</b> (from 38 to 97 cm)	SAA50402
<b>Abdominal belt kit standard size</b> (from 53 to 112 cm)	SAA50400
<b>Abdominal belt kit large size</b> (from 78 to 137 cm)	SAA50401
<b>Extension with abdominal belt</b> (recommended for the obese - L. 30 cm)	SAA50410
<b>Pelvic strap</b>	SAA50420

## Limb strap



The kit contains 2 limb strap, 1 magnetic key et 4 axis + magnetic bolt and twine. The limb strap are equipped with a velcro closure system allowing quick and easy control wrists and ankles in an emergency situation.

Description	Reference
<b>Limb strap kit small size</b>	SAA50301
<b>Limb strap kit standard version</b>	SAA50300