

▶ **Neop**

Medication and care
management system




PRATICIMA
healthcare solutions

▶ Neop

designed to cover all your needs

A range of trolleys and cabinets for the transport, storage and distribution of medicine and medical products, care delivery and for care unit equipment.



Design

- **Unique design** to meet all equipment needs of healthcare establishments
- **Asymmetric shape:** distinctive look
- **Harmonious colours,**
3 available colours:

Magenta pink

Apple green

Sky blue



Robustness

- **Monobloc box** made of aluminium or steel according to models: **perfect corrosion resistance**
- **Very high impact resistance,** outside bumper for a complete trolley protection
- **Qualitative castors,** to support heavy loads
- **+ Transfer range allowing strapping at all levels**



Ergonomics

Top shelf containing ergonomic gripping area:

perfect handling and workability

Ergonomic handle*:

natural position of hands and wrists and prevents shoulder strains

Opening / closing of roller shutters without bending down thanks to intermediate grip areas:

limits MSD**



Extension shelf*

to enlarge working surface

Motorisation*:

limits pushing efforts and prevents MSD**

*As an option

**MSD: Musculo-Skeletal Disorders



Hygiene

- **Monobloc box design:**
no retention zone
- **Resistant to disinfecting products of hospital environment**
- **100 % cleaning of the trolley:**
extractable roller shutter (without tool) and removable inside side panels




Customization

- **Wide range of internal and lateral accessories** compatible with full Neop range
- ➕ **Ward name identification** to keep furniture assignment in each ward
- ➕ **Codelock** (with spare key), **badge lock** or **badge lock with access control**


Range of products

- Distribution and care trolleys
- Emergency trolleys
- Mobile storage cabinets
- Transfer trolleys and cabinets
- Tubular or suspended working surface
- Mobile units
- Mobile worktops
- Wall units...



3 dimensions

600 x 400
400 x 400
400 x 600



6 heights

60, 80,
100, 120,
160, 180 cm



Find our product data sheets
on our website!

► PRATICIMA serving your projects

- **Project support:** sales force of 20 people specialized in each market segment (health and social care services)
- **Short delivery times**
- **Integrated logistics service,** delivery by our own trucks and drivers
- **Installation service** upon request
- ⊕ **Management process optimization of pharmaceutical products and consumable storage:** to increase the efficiency of tasks linked to internal supply management
- ⊕ **Asset audits and maintenance operations of modular material in healthcare facilities:** to guarantee an optimum service level and optimize investments



PRATICIMA
healthcare solutions

Rue du Pou du Ciel
ZI de Reyrieux – CS20129 Reyrieux
01600 TRÉVOUX – FRANCE
Tel: +33 (0)4 74 00 98 11
Fax: +33 (0)4 74 00 98 36
E-mail: praticima@praticima.fr

www.praticima.fr

