

Transfer Range

Rollbord Ecolite® transfer boards

The Rollbord Ecolite® transfer boards enable to transfer the patients without lifting them. Work conditions of medical staff can therefore be considerably improved.

As they are in accordance with hygiene standards, they are safe for skin and also very comfortable for the patients.

With their cover with silicone coating, The Ecolite® transfer boards are designed for a standard use.



Hygienic



Ergonomic

Advantages :

- ▶ Safe transfer of the patients
- ▶ Prevent medical staff from back injury
- ▶ Antibacterial, antifungal and fireproof coating
- ▶ Easy to clean and disinfect

Transfer range

Rollbord Ecolite® transfer boards

Technical specifications

- ▶ Polyethylene foam core
- ▶ Polyester sliding cover with Silicone coating
- ▶ Polyethylene ends
- ▶ Medical device class I according to 2017/745 European Regulation



Description	Use	Dimensions (in mm)	Weight	Capacity	Reference
« Ecolite® 90 » non foldable board	Operating theatre for urology and gynecology, emergencies, intensive care units...	L. 900 x 500	1,3 kg	150 kg	SAA71013
« Ecolite® 120 » foldable board	Radiology, operating theatres,...	Unfolded : 1200 x 400 Folded : 600 x 400	1,4 kg	150 kg	SAA71014
« Ecolite® 180 » foldable board	Radiology, operating theatres, intensive care units...	Unfolded : 1800 x 500 Folded : 870 x 500	2,8 kg	160 kg	SAA71012

Optional equipments

Description	Photo	Dimensions (in mm)	Reference
Wall-mounted support for transfer board made of wire stainless steel		L. 410 x D. 70 x H. 600	SAA40250
Wall-mounted roll support for single use covers made of wire stainless steel		L. 690 x D. 210 x H. 130	SAA40200
Case of 6 rolls of 50 single use covers		L. 2000 x D. 650	SAA40300