

Transfer range

# Rollbord Hightec® Transfer Boards

The Rollbord Hightec® transfer boards enable to transfer the patients without lifting them. Work conditions of medical staff can therefore be considerably improved.

As they are in accordance with hygiene standards, they are safe for skin and also very comfortable for the patients.

With their wear-resistant cover, the Rollbord Hightec® transfer boards are designed for intensive use.



Hygienic



Ergonomic

## Advantages :

- ▶ Safe transfer of the patients
- ▶ Prevent medical staff from back injury
- ▶ Antibacterial, antifungal and fireproof coating
- ▶ Easy to clean and disinfect
- ▶ Can cover spaces until 20 cm

Transfer range

# Rollbord Hightec® Transfer Boards

## Technical specifications

- ▶ Polyethylene foam core
- ▶ Polyester sliding cover with Teflon coating
- ▶ ABS handles integrated in the ends
- ▶ Medical device class I according to 2017/745 European Regulation



Description	Use	Dimensions (in mm)	Weight	Capacity	Reference
"Mini-gyn" non foldable board	Operating theatre for urology and gynecology, emergencies, intensive care units...	L. 910 x D. 500	1,4 kg	180 kg	SAA71005
"Standard" foldable board	Radiology, operating theatres, intensive care units...	Unfolded : L. 1770 x D. 500 Folded : L. 870 x D. 500	2,7 kg	240 kg	SAA71001
"Icu-bariatric" non foldable reinforced board	Intensive care units, transfer of bariatric people	L. 1770 x D. 500	3,8 kg	310 kg	SAA71004

## Optional equipments

Description	Photo	Dimensions (in mm)	Reference
Wall-mounted support for transfer board made of wire stainless steel		L. 410 x D. 70 x H. 600	SAA40250
Wall-mounted roll support for single use covers made of wire stainless steel		L. 690 x D. 210 x H. 130	SAA40200
Case of 6 rolls of 50 single use covers		L. 2000 x D. 650	SAA40300