


## Transfer range

# Transfer boards with curved design

The transfer boards with curved design are a good help for the transfer of a patient and are especially designed for people who have difficulties to stand up. They offer assistance and enable to change from one position to another without any efforts : transfer from bed to wheelchair, from chair to toilets or from chair to car.

Their « banana » curved design has been studied in cooperation with therapists in order to facilitate user's position and assure a stable and safe support.

These boards, which have already been sold in several tens of thousands units, are ideal for medicine wards, old people's homes, and home care.



Hygienic



Ergonomic

## Advantages

- ▶ Curved design (banana)
- ▶ Non-slip adhesive to avoid gliding
- ▶ Robust and light
- ▶ Can support until 150 kg
- ▶ Easy cleaning

## Transfer range

# Transfer boards with curved design

### « Banana » Transfer board

- ▶ Dimensions : L. 680 x D. 230 mm
- ▶ Weight : 1,2 kg
- ▶ Capacity (for a transfer with a distance lower than 20 cm) : 150 kg
- ▶ Material : HDPE resistant to acids, non toxic
- ▶ Easy to clean and to disinfect
- ▶ A (blue) non-slip adhesive guarantees the positioning on bed, wheelchair...

#### Description

#### Référence

« Banana » Transfer board

SAA72003



### « The biggest » Transfer board

- ▶ Dimensions : L. 1000 x D. 340 mm
- ▶ Weight : 2,9 kg
- ▶ Capacity (for a transfer with a distance lower than 30 cm) : 130 kg
- ▶ Material : HDPE resistant to acids, non toxic
- ▶ Easy to clean and to disinfect
- ▶ A (blue) non-slip adhesive guarantees the positioning on bed, wheelchair...

#### Description

#### Référence

« The biggest » Transfer board

SAA72001

